

## Axiom<sup>®</sup> Salmon Genotyping Array

The world's first-of-its-kind expert-designed Atlantic salmon (*Salmo salar*) genotyping array for high-resolution genetics research and aquaculture breeding

Axiom<sup>®</sup> Salmon Genotyping Array for Atlantic salmon (*Salmo salar*) was designed using Affymetrix<sup>®</sup> Expert Design Program in collaboration with The Roslin Institute (University of Edinburgh), Edinburgh Genomics, Landcatch Natural Selection Ltd., University of Stirling, and University of Glasgow in the United Kingdom.

### Benefits of high-density genotyping of salmon

Axiom Salmon Genotyping Array provides extensive information on polymorphisms across the salmon genome. Such dense SNP information is required to capture population-wide linkage disequilibrium for effective salmon GWAS.<sup>1</sup> Further, dense marker maps are required for effective development of breeding programs that will have sufficient accuracy for genomic selection.<sup>2</sup>

The high marker density on Axiom Salmon Genotyping Array ensures a broad coverage of the salmon genome and avoids underrepresentation of important regions of the salmon genome. This is a significant improvement over the existing Atlantic salmon genotyping tools.

Salmonids experienced a recent whole-genome duplication event, which caused Atlantic salmon to have a mosaic genome. Although most of the genome segregates in a diploid fashion, it can be difficult to distinguish genuine segregating SNPs from paralogous variation. Axiom Salmon Genotyping Array and corresponding Axiom<sup>®</sup> analysis tools partially overcome the challenges arising from this complexity, including false variant calls, cluster compression, and species-specific polymorphic site and structural diversity.

### Highlights

#### Diverse content:

>130,000 SNPs evenly spaced across the genome

- Preference given to transcribed regions
- Polymorphisms discovered using three different sequencing methodologies
- 4,714 markers from public sources
- 87 Y chromosome-specific probes for accurate gender determination

#### Multiple breeds represented:

- Farmed Scottish and Norwegian and wild Scottish, Irish, Norwegian, and Spanish populations

#### Demonstrated use:

- Detection of population structure
- Determining gender of salmon
- Genome-wide association study (GWAS) for resistance to sea lice<sup>3</sup>

#### Applications

##### Complex trait research:

- GWAS
- Fine mapping of QTL affecting economically important traits in salmonids

##### Molecular breeding:

- Association mapping and genomic selection
- Improved accuracy in aquaculture breeding programs
- Gender determination via Y chromosome-specific probes

##### Population and evolutionary genetics:

- Distinguish between fish of different origins
- Differentiate farmed and wild populations

##### SNP discovery and marker selection

SNP discovery was facilitated by sequencing both farmed and wild fish. Sequences were aligned to the draft Atlantic salmon reference genome (NCBI Assembly GCA\_000233375.1).

Due to genome duplication and mosaicism in the salmon genome, paralogous sequence variants (PSVs)<sup>4</sup> may interfere with allelic SNPs and complicate SNP discovery. Therefore, three distinct genome complexity reduction approaches were applied: reduced representation (RR-Seq), restriction site-associated DNA (RAD-Seq), and mRNA (RNA-Seq). The high sequence coverage and use of both pedigreed samples and a haploid individual were adopted in order to identify and remove PSVs and include informative polymorphisms. See Table 1 for details of these experiments.

**Table 1:** Techniques of DNA and RNA sequencing experiments used for SNP discovery

	RR-Seq	RAD-Seq	RNA-Seq
<b>Samples</b>	*Farmed (40), wild (16), haploid (1)	*Farmed (160)	*Farmed (72)
<b>Initial discovered SNPs</b>	472,072	467,268	816,570
<b>SNPs selected for array**</b>	99,097	83,151	229,754
<b>Final SNPs for array design***</b>	73,800	54,197	156,979

\*Farmed fish were representative samples of the Atlantic salmon breeding company Landcatch Natural Selection Ltd.

\*\*SNPs selected after filtering for PSVs, MAF, and repeat elements

\*\*\*SNPs selected after calculation of *in silico* design scores

### RR-Seq:

- One haploid individual was included to identify heterozygous PSVs, which were then excluded from the final array design
- SNPs with minor allele frequency (MAF)  $\leq 0.1$  were removed
- SNPs with known repeat elements were removed

### RAD-Seq:

- SNPs with >1 Mendelian errors across families were removed
- Repeat-masked SNPs were removed

### RNA-Seq:

- Repeat-masked SNPs were removed

### References

- Houston R., *et al.* Development and validation of a high-density SNP genotyping array for Atlantic salmon (*Salmo salar*). *PAG XXII*.
- Meuwissen T. H., Hayes B. J., Goddard M. E. Prediction of Total Genetic Value Using Genome-Wide Dense Marker Maps. *Genetics* **157**(4):1819–1829 (2001).
- International Aquafeed* **16**(1):3 (2013). [http://issuu.com/international\\_aquafeed/docs/iaf1304\\_w1](http://issuu.com/international_aquafeed/docs/iaf1304_w1)
- Fredman, D., *et al.* Complex SNP-related sequence variation in segmental genome duplications. *Nature Genetics* **36**(8):861–866 (2004).
- Yano A., *et al.* An immune-related gene evolved into the master sex-determining gene in rainbow trout, *Oncorhynchus mykiss*. *Current Biology* **22**(15):1423–1428 (2012).

### Ordering information

Part number	Product name	Description
550540	Axiom® Salmon Genotyping Array (ssalar01)	Contains one plate with 96 arrays. Reagents and GeneTitan® MC consumables must be quoted separately.
901606	Axiom® GeneTitan® Consumables Kit	Contains all GeneTitan® Instrument consumables required to process one array plate
901758	Axiom® 2.0 Reagent Kit	Includes all reagents (except isopropanol) for processing 96 DNA samples

Affymetrix, Inc. Tel: +1-888-362-2447 ■ Affymetrix UK Ltd. Tel: +44-(0)-1628-552550 ■ Affymetrix Japan K.K. Tel: +81-(0)3-6430-4020

Panomics Solutions Tel: +1-877-726-6642 [panomics.affymetrix.com](http://panomics.affymetrix.com) ■ USB Products Tel: +1-800-321-9322 [usb.affymetrix.com](http://usb.affymetrix.com)

[www.affymetrix.com](http://www.affymetrix.com) Please visit our website for international distributor contact information.

**For Research Use Only. Not for use in diagnostic procedures.**

P/N GGNO3229 Rev. 2

©Affymetrix, Inc. All rights reserved. Affymetrix®, Axiom®, Command Console®, CytoScan®, DMET™, GeneAtlas®, GeneChip®, GeneChip-compatible™, GeneTitan®, Genotyping Console™, myDesign™, NetAffx®, OncoScan™, Powered by Affymetrix™, PrimeView™, Procarta®, and QuantiGene® are trademarks or registered trademarks of Affymetrix, Inc. All other trademarks are the property of their respective owners.

Products may be covered by one or more of the following patents: U.S. Patent Nos. 5,445,934; 5,744,305; 5,945,334; 6,140,044; 6,399,365; 6,420,169; 6,551,817; 6,733,977; 7,629,164; 7,790,389 and D430,024 and other U.S. or foreign patents. Products are manufactured and sold under license from OGT under 5,700,637 and 6,054,270.

Putative SNPs were submitted to Affymetrix to calculate *in silico* design scores, which predict the likelihood of success in the Axiom® assay. SNPs with high scores were included for the final design.

Y chromosome-specific SNPs were included by conducting a homology search of the sdY gene, which was recently discovered in rainbow trout.<sup>5</sup> Partial sequences of the gene were validated with PCR, and two contigs were produced. Repetitive elements were removed, and 87 probes were designed for the array.

### Results

Axiom® Salmon Genotyping Array has been successfully used to genotype salmon samples of diverse origin from three main categories: farmed Scottish, farmed Norwegian, and wild fish from Scotland, Ireland, Norway, and Spain.

The Axiom Salmon Array has been used in studies that confirmed the sex locus in the Atlantic salmon with 100% concordance and showed that the array can be used to distinguish between different populations of both farmed and wild salmon. It is currently being applied in a GWAS for resistance to sea lice.

Samples were assayed on Axiom Salmon Genotyping Array using the standard Axiom® 2.0 assay and GeneTitan® Multi-Channel (MC) Instrument running GeneChip® Command Console® Software (AGCC) 3.2.4 or higher.

The data analysis and clustering were automated using Affymetrix' Genotyping Console™ Software 4.1.4 or higher and Affymetrix' SNPolisher™ package 1.4. Sample QC was performed as per the *Best Practice Supplement to Axiom® Genotyping Solution Data Analysis* (P/N 703083). SNP QC was performed as per the *SNPolisher™ User Guide*.