





Figure 2: PorcineSNP60 v2 BeadChip Probe Spacing—The PorcineSNP60 v2 BeadChip provides uniform coverage across the porcine genome.

### High-Quality Data

The 64,232 SNPs on the PorcineSNP60 v2 BeadChip were subjected to rigorous functional testing on multiple breeds to ensure strong performance using the Infinium HD Assay. High call rates and accurate genotype calls are required for successful whole-genome association studies. Because complex traits often have relatively small gene effects, potential associations can be missed if the assayed SNP in linkage disequilibrium with the SNP of interest has a low call rate or incorrect genotype call. Illumina ensures that every PorcineSNP60 v2 BeadChip offers > 99% average call rate across common breeds.

Table 1 shows the results from internal validation testing of the PorcineSNP60 v2 BeadChip content using samples provided by collaborators and the Porcine HapMap Consortium. Illumina scientists and collaborators identified and retained 2596 loci that appear to have an adjacent polymorphism or deletion among the breeds sampled. Although these loci yielded lower call rates when compared to most loci on the panel, they were retained because they might provide biologically relevant information, especially between breed groups. These performance and content validation results clearly demonstrate the high data quality delivered by the product. With such high data quality, the PorcineSNP60 v2 BeadChip provides researchers the highest accuracy and reliability for interrogating porcine genotypes in numerous porcine breeds.

Table 1: Performance and Specifications

Parameter	Results	Product Specification
Average call rate <sup>a</sup>	99.8%	> 99%
Reproducibility	100%	> 99.9%
Mendelian inconsistencies	0.08%	< 0.1%

a. Based on genotyping 554 reference samples.

### Illumina Solutions for Genotyping

The PorcineSNP60 v2 BeadChip is compatible with the iScan® and HiScan® Systems. These array scanners feature high-performance lasers and powerful optical systems that enable rapid scan times and precise assay detection.

The convenient modular design enables researchers to build out the system easily for evolving research needs. An optional Laboratory Information Management System (LIMS) is available to track samples accurately and efficiently. Robotic automation capabilities can be added to improve throughput for labs processing large numbers of samples. With the Infinium Assay workflow, data are processed directly into Illumina GenomeStudio® software to provide streamlined genotype calling, analysis, and reporting. Researchers can also choose to use the convenient FastTrack™ Genotyping service to have samples genotyped and data delivered in a format suitable for GWAS or QTL analysis.

### Summary

Developed through collaboration between Illumina scientists and leading thought leaders, the PorcineSNP60 v2 BeadChip features 64,232 evenly spaced SNPs that provide comprehensive coverage of the porcine genome, enabling a diverse range of genomic research applications. This 24-sample BeadChip, along with the proven Infinium HD assay, presents a powerful high-throughput solution for whole-genome studies for many porcine breeds.



